

ReWalk Discharge Plan

What Makes a Good Companion?

- Lives with or near the User
- Available on a regular basis both during and after the training (1hr, 1-3X a week)
- Physically able to provide assistance as needed (donning, regaining balance)
- Physically able to move (not necessarily lift) 60lbs
- Ability to communicate openly and calmly with the User

Potential Companions: (Only one is required, more are encouraged)

<u>Name</u> _____	<u>Relationship</u> _____	<u>Possible Issues to Overcome</u>
_____	_____	_____
_____	_____	_____
_____	_____	_____

Transport Plan: Dolly _____

Car _____

Other _____

Storage Plan: Bench _____

Charging _____

Other _____

Usage Plan: Frequency/Duration _____

Location _____

Checklist for ReWalk Device Home Release:

Before the User takes the device home, please complete the following:

- Home visit by a training PT is ideal, if not possible, discuss home environment and ReWalk use at home and in the community
- Sign and collect signatures for Basic and Advanced inventories- email a copy to your ReWalk contact person
- Provide Field Service Engineer contact information
- Cover the charging requirements
 - Demonstrate how to open the waistpack to disconnect batteries from INF
 - P6 Quick User Reference has charging instructions
- Demonstrate how to use the computer and connect to ReWalk program
- Demonstrate tightening of bolts with a torque wrench
- Ankle joint instructions:
 - Ask the user to keep an eye on ankle joints and instruct to use allen wrench to make the adjustment if needed.
- Provide a printed copy of P6 quick user reference (DOC0621-14)
[Link to the document](#)
- Provide information about a dolly and a bench.

Example of a dolly:

[Amazon link for dolly](#)

Example of a foldable bench:

[Amazon link to a bench](#)